



Gamma Spectroscopy of Natural Radionuclides: techniques and instrumentation

Yuri Venturini

June 24th, 2025

Carpathian Summer School, Sinaia - RO



Link to CAEN Educational



CAEN Edu introduction


CAEN S.p.A. is an industrial spin-off of the Italian Institute of Nuclear Physics founded in 1979

Core business: Electronic Instrumentation for High-energy physics

→ High Voltage and Digital Pulse Processing

- **HQ** in Italy
- **3** direct sales offices in U.S., Germany, India
- **30+** distributors worldwide
- **160+** specialized employees
- **5000+** customers

- *High Energy Physics*
- *Astrophysics*
- *Neutrino Physics*
- *Dark Matter Investigation*
- *Nuclear Physics*
- *Material Science*
- *Medical Applications*
- *Homeland Security*



Modern physics experiments for Advanced Labs
based on the latest technologies and methods




80+
Activities



14
Edu kits



50+
Institutions



Mission

- ❖ Guide students towards the comprehension of different physics phenomena with experiments based on **state-of-the art technologies**, instruments and methods
- ❖ Target the experiment depending on the **student educational level**. From high school level (grade 11,12) up to undergraduate physics laboratory and PhD courses.

Nuclear Physics and Radioactivity

γ Spectroscopy

- ✓ Detecting γ -Radiation
- ✓ Poisson and Gaussian Distributions
- ✓ Energy Resolution
- ✓ System Calibration: Linearity and Resolution
- ✓ A comparison of different scintillating crystals: Light Yield, Decay Time and resolution
- ✓ γ -Radiation Absorption
- ✓ Photonuclear cross-section/Compton Scattering cross-section

β -Radiation

- ✓ Response of a Plastic Scintillating Tile
- ✓ β Spectroscopy
- ✓ β -radiation: Transmission through Matter
- ✓ β -Radiation as a Method to Measure Paper Sheet and thin layer thickness

Nuclear Imaging - PET

- ✓ Basic Measurements: γ Spectroscopy and System
- ✓ Positron Annihilation Detection
- ✓ Two-dimensional Reconstruction of Source
- ✓ Spatial Resolution

γ Environmental Radioactivity (outdoor)

- ✓ Environmental monitoring in land field
- ✓ Ground Coverage Effect on the Environmental
- ✓ Human Body Radioactivity
- ✓ Environmental detection as a function of the s
- ✓ Radioactivity maps production
- ✓ Radiological evaluation of the building materia
- ✓ Geochemical and mineral exploration

γ Environmental Radioactivity (indoor)

- ✓ Energy calibration of System based on LYSO crystal
- ✓ Background Measurements
- ✓ Fertilizer and photopeak identification
- ✓ Identifications Sample Test
- ✓ Soil sample identification
- ✓ Samples Comparison
- ✓ Radon passive measurements



GM Detectors

- ✓ Statistics: Uncertainty as a function of live time
- ✓ Environmental Background
- ✓ Lead Shielding Effect on Environmental Radioactivity
- ✓ Detecting Ionizing-Radiation
- ✓ Samples Comparison



Particle Physics

Photons

- ✓ Quantum Nature of Light
- ✓ Hands-on Photon Counting Statistics



Cosmic Rays

- ✓ Statistics
- ✓ Muons Detection
- ✓ Muons Spectrum
- ✓ Muons Vertical Flux on Horizontal Detector
- ✓ Zenith Dependence of Muons Flux
- ✓ Random Coincidence
- ✓ Detection Efficiency
- ✓ Cosmic Flux as a function of the altitude
- ✓ Cosmic Shower Detection
- ✓ Environmental and Cosmic Radiation
- ✓ Absorption Measurements
- ✓ Solar Activity Monitoring



Particle Detector Characterization

Silicon Photomultiplier (SiPM)

- ✓ SiPM Characterization
- ✓ Dependence of the SiPM Properties on the bias voltage
- ✓ Temperature Effects on SiPM Properties




Photomultiplier Tube (PMT)

- ✓ Measurement of Photomultiplier Plateau Curves



Pulse Processing: Open FPGA

- ✓ Analog signal acquisition and waveform Visualization
- ✓ Waveform digitizer with leading edge trigger.....



Short Guide

Main Topics:

- Experiment task
- Short description
- Equipment list
- Requirements
- Quick guide
- Experimental results



Detailed Guide

Guide Topics:

- General Information
- Introduction
- Physics Pills
- Required Equipment
- Getting Started
- Experimental Procedure
- Results
- Links related to this topic

Theory
(just 10 minutes 😊)

Commonly the concept of **Radioactivity** is always associated with a **dangerous feeling**!



Nuclear Wars



Nuclear accidents




Radioactive Waste

Radiation is a **natural part** of our environment!

Radioactive sources:

- **Natural: NORM** (Naturally Occurring Radioactive Material), soil, water, air and food contribute to our exposure to ionizing radiation
- **Cosmic rays**
- **Industrial:** nuclear elements produced by industry
- **Medical:** nuclear medicine



https://www.who.int/ionizing_radiation/env/en/

Radioisotope types

Natural radionuclides [mean abundance]:

- ^{40}K [2-2.5] %
- $^{238,235}\text{U}$ [2-2.5] ppm
- ^{232}Th [8 -12] ppm
- All has an half life $T_{1/2}$ > than the age of the solar system

Cosmic rays

Muons -> The intensity depends on the altitude and direction (N/S vs E/W)

Radionuclides from cosmic rays

- $^{14}\text{C}, ^7\text{Be}, ^3\text{He}$

Energy Distribution in Cosmic Rays
W. G. Pollard
Phys. Rev. **44**, 703 – Published 1 November 1933

Artificial radionuclides

- From bombs or Nuclear power plants (ex. ^{137}Cs , actinides)
- Industrial ($^{133\text{m}}\text{Xe}$, ^{133}Ba , ^{241}Am) and medical (^{19}F , ^{67}Ga) radioisotopes

During the creation of the Earth, most of the elements initially produced were radioactive and they have been decayed to more stable forms.

The original radioactive elements still present on Earth are those that have a halftime comparable to the Earth. They are responsible for environmental radioactivity and internal warming of the planet and originate from elements very heavy without stable isotopes.

They mostly decay through the α and β channels

Element	Radioisotopes	Isotopic Adundance	Half time	Tipycal Adundance
Potassium	^{40}K	0.012%	1.3×10^9 years	0.02 g/g [2%]
Uranium	^{238}U	99.3 %	4.5×10^9 years	2-3 $\mu\text{g/g}$ [ppm]
Thorium	^{232}Th	100 %	14.1×10^9 years	8-12 $\mu\text{g/g}$ [ppm]

The Periodic Table of the Elements

Notes:

- as of yet, elements 113-118 have no official name designated by the IUPAC.
- 1 kJ/mol = 96.485 eV
- all elements are implied to have an oxidation state of zero.




Most of the matter nuclei are stable even when they undergo chemical reactions

Some nuclei, however, are unstable (**radionuclides**): they transform spontaneously and **reach stability by emitting radiation**.

This process is called **radioactive decay**!

Nuclear decay is a random event... like popcorn popping!







X = Element symbol Z = Atomic Number (number of protons) N = Neutronic number (number of neutrons)

A = Mass Number (Z + N)

Isotopes

Elements with the same atomic number and different mass number

Ex: Hydrogen isotopes



The **activity** is defined as the number of decay per second. The international system unit is the Becquerel [Bq]

It is an exponential law:


$$N(t) = N_0 e^{-\lambda t}$$

Where:



- λ is the **decay constant** of the nuclide
- **Mean time** of the element: $\tau = 1/\lambda$
- **Half life**: $t_{1/2} = \tau \lg 2 \approx 0,7\tau$

The radioactive decay is a physics phenomena happening when an instable nucleus reaches a new state of equilibrium emitting particle or radiation


Type	Mass	Q	Description
α	$6.68 \cdot 10^{-27}$ Kg	+2	${}^4_2\text{He}$
β	$0.9 \cdot 10^{-30}$ Kg	-1	Electron (e^-)
γ	0	0	Electromagnetic radiation (photon)



- In the **α and β decays**, the nucleus is transformed into a different type emitting alpha or beta particles




- In the **γ decay**, the nuclide emits photons to go to a lower energy and more stable state. The photon energy corresponds to the excess of energy between the two states





Usually, the gamma decay is generated after an alpha or a beta decay which generate excited states. The **excited states** goes to a more stable state emitting photons and without changing state

The radioactive decay is a physics phenomena happening when an instable nucleus reaches a new state of equilibrium emitting particle or radiation

Type	Mass	Q	Description
α	$6.68 \cdot 10^{-27}$ Kg	+2	${}^4_2\text{He}$
β	$0.9 \cdot 10^{-30}$ Kg	-1	Electron (e^-)
γ	0	0	Electromagnetic radiation (photon)



- In the **α and β decays**, the nucleus is transformed into a different type emitting alpha or beta particles




- In the **γ decay**, the nuclide emits photons to go to a lower energy and more stable state. The photon energy corresponds to the excess of energy between the two states

Natural radioisotopes are primarily gamma emitters



Usually, **γ decay** is generated after an alpha or a beta decay which generate excited states. The **excited states** goes to a more stable state emitting photons and without changing state




Potassium is essential for living. In the human body, most of the potassium is stored in the muscles. Potassium is also present in the soil, building materials, plants, animals and it is used in fertilizers. In nature exists 3 potassium isotopes:

Type	^{39}K	^{40}K	^{41}K
n° protons	19	19	19
n° neutrons	20	21 	19
Abundance (%)	93.26	0.01	6.73
$t_{1/2}$	Stable	$1.3 \cdot 10^9$ y	stable


Considering a hundred thousand potassium atoms, only 12 are actually radioactive!



^{40}Ar decays emitting a gamma of 1460 keV (10%)
Measurement % \rightarrow 1% = 313 Bq/kg



The ^{238}U is the most common isotopes of the uranium element with a relative abundance of 99% and an half life of about 4,5 billions of years. The decay chain of an element down to its fundamental state is also called **radioactive series**.





- All the decay products found in the chain have a shorter average life in comparison to the generating elements of the series.
- The **secular equilibrium** is the situation when the quantity of an isotope remains constant having the same production and decay rate. This equilibrium is usually broken when one of the sons is a gas that goes away.

14 transformations are needed for the ^{238}U to be transformed into the stable ^{206}Pb

^{232}Th has 142 neutrons, is the most stable isotopes of the thorium family (10^{10} years) and represents almost all thorium existing in nature.



- It decays slowly with alpha decay on ^{198}Ra . The thorium series finishes with a stable element: the ^{208}Pb . This series contains the following elements: astatine, bismuth, lead, polonium, radium and radon
- The quantity of thorium in the Earth is 3 or 4 times larger than uranium. It can be extracted from monazite sands or as a sub-product of rare Earth elements extraction



Instrumentation and measurements

Energy deposition by an ionizing particle:

- Generation of light
- Transmission of scintillation light
- Detection



What are scintillators used for?

- To measure the energy released
- To measure the passage time of radiation

Organic (plastics or liquid solutions)

Up to 10000 photons per MeV

Low Z

$\rho \sim 1 \text{ g/cm}^3$

Doped, choice of emission wavelength

ns decay times

Relatively inexpensive

Inorganic (crystalline structure)

Up to 40000 photons per MeV

High Z

Large variety of Z and ρ

Un-doped and doped

ns to μs decay times

Expensive



Thallium doped sodium iodide, **NaI(Tl)**, is the most widely used scintillation material, it has the greatest light output and convenient emission range

NaI(Tl)


General Properties

Density(g/cm ³)	3.67
Melting point(K)	924
Wavelength of emission peak(nm)	415
Light output(Photons/Mev)	40,000
Decay time(ns)	264
Cleavage plane	(100)
Hygroscopic	Yes
Refractive index	1.85
Hardness(Mho)	2



Photodetector → From photons to electric current!

Photomultiplier Tubes (PMT) are composed of a photocathode, collection optics and multiplier section. The overall electrical signal is collected at the anode.



❑ PMT collects and transforms the light produced by the scintillator into an electrical signal

❑ The intensity of the output current pulse is proportional to the energy of the incident photon!






CAEN Gamma *stream* [S2580] is a compact and portable system for gamma ray spectroscopy with scintillation detectors, which provides an active **Multi-Channel Analyzer** (MCA) integrated in a 14-pin photo-multiplier tube (PMT) base.

Gamma *stream* fully integrates in a stand-alone device the high voltage to bias the PMT, the preamplifier to shape the signal from detector, and the MCA for a complete Pulse Height Analysis online.

Gamma *stream* makes easy the measurements with scintillation detectors **NaI(Tl)** [0.3l] with no need of additional cables.

- High Voltage Power Supply ($0 \div +1500\text{V}/500 \mu\text{A}$)
- Charge Sensitive Preamplifier
- digital Multi-Channel Analyzer (12-bit and 62.5 MHz ADC) for scintillation spectroscopy
- Specialized for NaI(Tl), LaBr₃(Ce), and CeBr₃ with standard 14-pin and 10-8 stages PMTs
- Full stand-alone operation with embedded CPU, data storage (SSD) unit, and power supply for up to 6÷8 hours operation
- Wired and wireless connectivity via USB, Ethernet, Wifi and Bluetooth
- Acquisition modes: PHA, PHA with time stamp, Signal Inspector





The acquisition channel is proportional to the energy of the incident photons!

The MultiChannel Analyzer MCA classifies input pulses base on their height saving them in a memory and are associated to an ADC. The output of every channel can be visualized in a pulse amplitude spectra.

An Analog-to-Digital Converter (ADC) generates a digital signal proportional to the amplitude of an input pulse. Since these output pulses are proportional to the energies of the incident radiation, the ADC can be used combined to a MultiChannel Analyzer (MCA) to generate energy distributions (spectra) of radioactive samples.



Ex. Monochromatic photons beam



The **MultiChannel Analyzer MCA** classifies input pulses base on their height saving them in a memory and are associated to an ADC. The output of every channel can be visualized in a pulse amplitude spectra.



An Analog-to-Digital Converter (ADC) generates a digital signal proportional to the amplitude of an input pulse. Since these output pulses are proportional to the energies of the incident radiation, the ADC can be used combined to a MultiChannel Analyzer (MCA) to generate energy distributions (spectra) of radioactive samples.

Ex. Monochromatic photons beam




The **MultiChannel Analyzer MCA** classifies input pulses base on their height saving them in a memory and are associated to an ADC. The output of every channel can be visualized in a pulse amplitude spectra.


An Analog-to-Digital Converter (ADC) generates a digital signal proportional to the amplitude of an input pulse. Since these output pulses are proportional to the energies of the incident radiation, the ADC can be used combined to a MultiChannel Analyzer (MCA) to generate energy distributions (spectra) of radioactive samples.



The photopeaks characterize the gamma spectrum. Each photopeak corresponds to the photons collected by the detector with an energy value equal to the emission ones. These photons release all their energy into detector.



The photopeaks characterize the gamma spectrum. Each photopeak corresponds to the photons coming into detector with an energy value equal to the emission ones. These photons release all their energy into detector.



In the energy range of the environmental measurements the calibration in energy corresponds to a linear transformation

=> Knowing the energy of 2 peaks it is possible to extract the equation of the line from 2 points

$^{40}\text{K} \Rightarrow E_K = 1460 \text{ keV}$
 $^{208}\text{Tl} \Rightarrow E_{Th} = 2614 \text{ keV}$

^{208}Tl is coming from the radioactive chain of the ^{232}Th and is the highest energy gamma from natural sources


$$A = (E_{Th}, Ch_{Th})$$

$$B = (E_K, Ch_K)$$

$$\frac{Ch - Ch_K}{Ch_{Th} - Ch_K} = \frac{E - E_K}{E_{Th} - E_K}$$

$$\frac{Ch - Ch_K}{Ch_{Th} - Ch_K} = \frac{E - 1460 \text{ keV}}{(2614 - 1460) \text{ keV}}$$



Multichannel Calibration



Definition of the **Region(s) of Interest (ROI)** of the energy spectra. These windows are used to define the photopeak regions required to calculate the correspondant areas (integral of the ROI).

NOTE: In every ROI there are different contribution effects (photopeak, compton continuum, background continuum, etc...)!


Window	Isotope	Photopeak Energy (keV)	ROI (keV)
Potassio	^{40}K	1460	1370-1570
Uranio	^{214}Bi	1765	1660-1860
Torio	^{208}Tl	2614	2410-2810



Count rate: number of counts per unit of time

$$n_i[\text{cps}] = \frac{N_i[\text{conteggi}]}{T[\text{s}]}$$

Calibration: from counts to the abundancies



$$N(\text{cps}) = A \cdot S$$

$$S = A^{-1} \cdot N$$

S is the sensitivity coefficient of the detector


$$A = N \cdot S^{-1}$$

Calibration site characterized by known abundancies (A) of U, Th and U.

Hands-on
(the funny part 😊)




Google
Earth
Software






- Press the power button
- Verify that the status light is green
- Place the backpack at the point of interest





- Activate the tablet and connect it via Bluetooth to the instrument
- Connect the tablet to a WIFI network
- Launch GammaEDU App
- Select «New Acquisition»








Set the measurement name

Take note of the surrounding environment

Take a measurement situ picture

Set the acquisition time





CAEN, Viareggio

Date and time: 2022_05_02_17_25
Environment: parck pineta


Total activity concentration and Isotopic abundances			
Total Activity Bq/kg	K %	eU ppm	eTh ppm
525 ± 31	1,5 ± 0,1	1,5 ± 0,4	5,8 ± 0,9

Typical abundancies

^{40}K	0.02 g/g [2%]
^{238}U	3 $\mu\text{g/g}$ [ppm]
^{232}Th	10 $\mu\text{g/g}$ [ppm]

Gamma spectroscopy in situ

- distribution of natural or artificial radioisotopes
- study possible radiological contamination
- studies in the field of earth sciences



- The main contributor on the overall natural indoor radioactivity is **$^{222}/^{220}\text{Rn} \rightarrow \text{U-238, Th-232}$**
- Emissions can be measured through **photopeaks**:
 - K-40 \rightarrow **1460 keV**
 - Tl-208 (from Th-232) \rightarrow **2614 keV**
 - Bi-214 (from U-238) \rightarrow 609 keV, **1764 keV**
- Compare indoor radiation** to typical environmental background or soil values.

Student Tasks:



- Measure background vs wall/floor material spectra (5 minutes acquisition)
- Identify peaks \rightarrow Assign isotopes.
- Estimate activity by comparing counts to a reference sample (e.g., earth's crust)

	Isotopic Abundances		
	^{238}eU [ppm]	^{232}Th [ppm]	^{40}K [%]
<i>Reference Values Range</i>	[2 ; 2.5]]	[8 ; 12]]	[2 ; 2.5]]
<i>Tuff Dwelling (4° floor)</i>	10 ± 1	31 ± 1	6.9 ± 0.2
<i>Modern Building (1° floor)</i>	2.8 ± 0.6	8.8 ± 1.1	1.6 ± 0.1
<i>Country House (0° floor)</i>	6.8 ± 0.9	17.6 ± 1.6	3.4 ± 0.2

- Soil emissions are **attenuated by distance and air** → lower intensity as you move away.
- **Attenuation in air** (for extended sources like soil) → $I(d) = I(0) \cdot e^{-\mu d}$
 - $I(0)$ is the count rate at the surface
 - μ is the **linear attenuation coefficient** (gamma energy and medium)
- The **height of the detector** affects the **lateral horizon** – i.e., **how far around** the detector the signal originates. In practice, you observe an **initial steep drop**, then a flattening due to wide angular coverage of the detector and scattering.

Student Tasks:

- Acquire spectra at several heights: 0 cm, 0.5 m, 1.5 m
- Plot count rate vs distance.
- Fit the data





Thanks for your attention!

Yuri Venturini on behalf of CAEN SpA

SP5630EN
Environmental kit



SP5630ENP
Environmental kit Plus



Two solutions for
environmental
radioactivity
indoor!



Environmental
radioactivity
Outdoor!



SP5640
Gamma EDU



SP5640 – GammaEDU

I.C.
Don Lorenzo
Milani
Viareggio, Italy
May 2023



RadioLAB
IRSOIL&WATER
2023
Spring School
Calabria, Italy
April 2023




*Portable detection backpack for
environmental radioactivity!*

RadioLAB Project
Sicily, Italy
April 2023




**High School
Students**
Tuscany, Italy
May 2022







			Equipment									
Section	Subsection	Experiment	SP5600C	SP5600D	SP5600E	SP5600AN	SP5600EMU	SP5700	SP5701	SP5620CH	SP5630EN	★ SP5640
Nuclear Physics and Radioactivity	γ Environmental Radioactivity (outdoor)	Environmental monitoring in land field										★
		γ Environmental Detection as a function of the soil distance										★
		Radioactivity maps production									★	★
		Mapping of potential radon-prone areas									★	★
		Radiological evaluation of the building materials									★	★
		Geochemical and mineral exploration									★	★




Radiological evaluation of the building materials




Geochemical and mineral exploration



Radioactivity maps production




Human body Radioactivity




γ Environmental detection as a function of the soil distance




Environmental monitoring in field




Ground coverage Effect on the Environmental Monitoring




Soil water content evaluation with gamma ray spectroscopy




Photoelectric effect



Compton effect



Pair production



- Photon interacts ONLY with atomic electrons, No with a free electron
- Most probable electron from the outer electronic shells
- In the interaction, creation of an electron-ion couple. The hole is immediately filled with free electron capture or reorganization of the electronic structure => Following emission of x-rays or Auger electrons
- Typical effect of the low energy region in the ~ keV region

- Photon interacts atomic electrons of the absorber material
- The interacted photon is scattered with a certain angle (θ) and the energy transferred allows the electron emission from the atom
- All the emission angles are allowed, but the scattered angle is energy dependent => Klein-Nishina Formula
- Typical interaction of radioisotope energy, hundreds of keV region

- Threshold effect => The pair production is possible only if the energy of the photon is larger than twice the energy at rest of the electron
- It is possible only in the Coulombian field of the nucleus
- The photon disappear for appearing a couple made of electron and positron
- The energy in above the threshold energy goes into kinetic energy of the couple
- The positron is annihilated rapidly with the following emission of 2 photons
- Typical interaction in the several MeV region

SPARES

Nuclear Physics

- Statistics
- Environmental Radioactivity

From geology to daily life, **RockyRAD** bridges the gap: discover the fascinating world of rock radioactivity and then expand your horizon to detect the unseen radiation in our everyday surroundings.


Portable Geiger Counter for nuclear radioactivity radiation!

Each **RockyRAD kit** is enhanced with a set of rock samples from different origins, allowing students to immediately begin their detection experiments.



Android App under development!

- ✓ *Detector:* Geiger-Müller Tube
- ✓ *Display Information:* Total Counts, Counts Per Minute, Equivalent Dose Rate
- ✓ *Wi-Fi* for data download
- ✓ *Bluetooth* connection
- ✓ *Rechargeable Battery* (USB-C)



SP5630EN Environmental kit



S2570B i-Spector Digital
18x18mm² - ASSEMBLY
(CsI 18x18x30mm³)

Samples



Empty Beaker &
Test Sample



Fertilizer and
Rock Samples



Canisters of
Activated Carbon



Calibration Crystal
(Lu1.8Y.2SiO5:Ce)


SP5630ENP Environmental kit Plus



S2570D i-Spector Digital
18x18mm² - CsI ASSEMBLY



BGO scintillating
crystal



**Shielding
Kit**

Samples



Empty Beaker &
Test Sample



Fertilizer and
Rock Samples



Canisters of
Activated Carbon



Calibration Crystal
(Lu1.8Y.2SiO5:Ce)

SP5640 Gamma EDU



S2580 GAMMASTREAM
Digital MCA Tube Base



Nal(Tl) [0.31]
scintillation
detectors




Tablet 10"




SP5630EN – Environmental kit

*Environmental gamma radiation measurements
with SiPM based instrumentation!*



S2570-iSpector Digital

- The system is based on a SiPM area $18 \times 18 \text{ mm}^2$
All SiPMs of the area are connected in parallel to increase the active area of the matrix.
- It integrates a shaper, a peak stretcher and a peak ADC to implement a simple MCA (4K).
- Scintillator Crystal: CsI $18 \times 18 \times 30 \text{ mm}^3$
- Connectivity: Ethernet
- Software: Web GUI



Samples



Empty Beaker & Test Sample



Fertilizer and Rock Samples




Canisters of Activated Carbon



**Calibration Crystal
(Lu1.8Y2SiO5:Ce)**


LYSO (Lu_{1.8}Y_{0.2}SiO₅:Ce) (Cerium-doped Lutetium Yttrium Orthosilicate) Scintillating Crystal

Scintillator based on Lutetium (Lu) like **LSO** and **LYSO (Lu_{1.8}Y_{0.2}SiO₅:Ce)** are usually used in PET applications thanks to their high stopping power (high Z), high light yield and very short decay time (very fast signals). It is a non-hygroscopic scintillator.






2,6% of the natural Lu is ¹⁷⁶Lu, a radioisotope with a long half life decaying via two different beta decays followed by gamma emissions.

$$T_{1/2} \sim 3,6 \times 10^{10} \text{ y}$$



Properties	Value
Cleavage Planes	None
Decay Constant (ns)	40
Density (g cm ⁻³)	7.1
Emission Spectral Range (nm)	380-480
Melting Point (K)	2323
Peak Scintillation Wavelength (nm)	420
Photons/MeV	32000
Radiation Length (cm)	1.15
Refractive Index at Peak	1.81
Emission	
Solubility (g/100g H ₂ O @ 300K)	Insoluble
Stability	Good
Structure	Cubic




Test Sample

One of the radioactive objects of common use, in past especially.

Canisters of Activated Carbon for Radon Passive Measurements

- Diameter: 10 cm
- Height: 3 cm
- Content: 70 g of activated carbon



The activated carbons are enclosed in metal containers called "canisters".

Covered by a thin double-mesh metal mesh (diffusive barrier).





The diffusion barrier serves to eliminate the air flows inside the basket, which can favor the re-emission of radon.

The method consists in carrying out gamma spectrometry measurements on the baskets after the radon has been adsorbed by them!

After 6-7 days, the loss due to decay prevails over the accumulation by adsorption.


Features:

- ❖ useful for short-term measurements: 2 - 7 days
- ❖ strong dependence of the response on humidity




Gamma Energy lines:

- 295 keV and 352 keV from ^{214}Pb
- 609 keV from ^{214}Bi




- The systems are based on a SiPM area $18 \times 18 \text{ mm}^2$
Hamamatsu S14160-60520HS
- All SiPMs of the area are connected in parallel to increase the active area of the matrix.

Silicon Photomultiplier (SiPM) are high density (up to $10^4/\text{mm}^2$) matrix of diodes with a common output, working in Geiger-Müller regime. Each cell is a pixel with a binary signal. It can detect a single photon!



- High Gain
- Low Voltage
- High photon number resolving power
- Wide dynamic range
- Good timing capability
- Low cost
- Withstanding to magnetic field


The high uniformity of pixel structure guarantees no avalanche fluctuations



Linear response if the average number of photoelectrons/pixel is less than one


Number of pixel determines the SiPM **dynamic range**

Binary information



$$A = \sum A_i$$

Analog information






- The systems are based on a SiPM area $18 \times 18 \text{ mm}^2$
Hamamatsu S14160-60520HS
- All SIPMs of the area are connected in parallel to increase the active area of the matrix.
- They integrate a shaper, a peak stretcher and a peak ADC to implement a simple MCA (4K).
- Scintillator Crystal: CsI $18 \times 18 \times 30 \text{ mm}^3$
- Energy Range: 30 keV to 3 MeV
- Energy Resolution (FWHM): $<6 \%$ @ 662 keV
- Connectivity: Ethernet
- Software: Web GUI

TECHNICAL SPECIFICATIONS

Supply Voltage	8-13V (12 V typ.)
Power consumption	3W max.
Preamplifier bandwidth	$>1\text{GHz}$
Preamplifier gain	x 5
Shaping time	180 ns
Output signal	- 4 ... +4 V , 170 mA
HV Power supply	20-80 V (10mA)
HV accuracy	1 mV
Thermal feedback accuracy	$0.01^\circ\text{C} - 1\text{mV}$
MCA nr. of channels	4096



[i-Spector Digital block diagram](#)



- The systems are based on a SiPM area $18 \times 18 \text{ mm}^2$
Hamamatsu S14160-60520HS
- All SiPMs of the area are connected in parallel to increase the active area of the matrix.
- They integrate a shaper, a peak stretcher and a peak ADC to implement a simple MCA (4K).
- Scintillator Crystal: CsI $18 \times 18 \times 30 \text{ mm}^3$
- Energy Range: 30 keV to 3 MeV
- Energy Resolution (FWHM): $< 6 \% @ 662 \text{ keV}$
- Connectivity: Ethernet
- Software: Web GUI

TECHNICAL SPECIFICATIONS

Supply Voltage	8-13V (12 V typ.)
Power consumption	3W max.
Preamplifier bandwidth	$> 1 \text{ GHz}$
Preamplifier gain	x 5
Shaping time	180 ns
Output signal	- 4 ... +4 V , 170 mA
HV Power supply	20-80 V (10mA)
HV accuracy	1 mV
Thermal feedback accuracy	$0.01^\circ\text{C} - 1 \text{ mV}$
MCA nr. of channels	4096

Web-based GUI for unit control and data analysis

i-Spector Digital can be easily controlled through its dedicated web graphical user interface, with no needs to install a dedicated software. The user can configure the module and visualize the acquired spectrum.

Thanks to the internal circular memory buffer, the last 1-hour recording can then be downloaded by the web interface.






- The systems are based on a SiPM area $18 \times 18 \text{ mm}^2$
Hamamatsu S14160-60520HS
- All SIPMs of the area are connected in parallel to increase the active area of the matrix.
- They integrate a shaper, a peak stretcher and a peak ADC to implement a simple MCA (4K).
- Scintillator Crystal: CsI $18 \times 18 \times 30 \text{ mm}^3$
- Energy Range: 30 keV to 3 MeV
- Energy Resolution (FWHM): $<6 \%$ @ 662 keV
- Connectivity: Ethernet
- Software: Web GUI

CsI(Tl) information:

This scintillator offers a high light yield and emits at a wavelength very suitable for silicon photomultipliers (SiPMs). Typical applications include arrays of this material used in security imaging systems, such as baggage scanners.


BGO information:

A relatively hard, high density, non-hydroscopic crystal with good gamma ray absorption. Often used for PET imaging and high energy physics applications as Compton shields.




Properties	CsI	BGO
Cleavage Planes	None	None
Decay Constant (ns)	1000	300
Density (g cm^{-3})	4.51	7.13
Emission Spectral Range (nm)	350-725	350-650
Gamma and X-ray absorption coefficients (cm^{-1})	0.48 at 660keV 10.00 at 100KeV	-
Melting Point (K)	894	1323
Peak Scintillation Wavelength (nm)	550	480
Photons/MeV	52000	8500
Radiation Length (cm)	1.86	1.13
Refractive Index at Peak Emission	1.78	2.15
Solubility ($\text{g}/100\text{g H}_2\text{O}$ @ 300K)	44.0	Insoluble
Stability	Slightly Hygroscopic	Good
Structure	BCC	Cubic
Thermal Conductivity ($\text{W}\cdot\text{m}^{-1}\cdot\text{K}^{-1}$) @ 300K	1.13	-
Transmission Range (nm)	240-70000	470-7500


Hands-On...




1. **Unboxing and Assembling**
2. **Software Setup**
3. **Energy calibration of the system based on LYSO crystal (time base = 10')**
4. **Calibration verification and tuning with Potassium Chloride sample (time run = 30')**
5. Background measurement (time run = 30')
6. Rock sample Spectrum (time run = 30')
7. Test sample radiation identification
8. Analysis of spectra and superposition
9. Passive Radon Measurements

- **IP address** of the i-Spector for Ethernet connection is **192.168.50.2**
- Configure the Ethernet network of the PC from the “Network and Sharing Center”
- Open a web browser (**Microsoft Edge browser is suggested**) and enter the web address **192.168.50.2**. The homepage of the graphical web interface will open.







- **Trigger Threshold (LSB):** threshold for the derivative trigger
- **Post Trigger Inhibit (ns):** set the time after a trigger for which any other trigger is inhibited
- **Pre-Integration (ns):** set how much time before the trigger the charge integration is started
- **Integration (us):** set the charge integration gate
- **Gain:** set the energy digital gain to be applied to the spectrum
- **Pileup Inhibition (us):** set the time after the integration gate for which the acquisition of any other event acquisition is inhibited
- **Pileup Penalty (us):** set the trigger inhibition gate to be opened after a pile up
- **Baseline Inhibition (us):** set the time after the integration gate for which the baseline is not calculated
- **Baseline Length (samples):** set the number of samples used to calculate the baseline
- **Target Mode:** set the acquisition mode as Free Running or with a target in Time (ms) or Counts
- **Target Value:** set the target value in time or counts, accordingly to the Target Mode



- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO
- 3) Background measurement (time run = 30')



Scintillator




Energy deposition by an ionizing particle:

- Generation of light
- Transmission of scintillation light
- Detection

CsI(Tl) has a light output of 54 photons/keV and average decay time of about 1μs for γ-rays

Silicon Photomultiplier (SiPM)


Silicon Photomultiplier (SiPM) is detector made of a matrix of silicon cells (diodes). Each diode is a pixel with a binary signal. It can detect a single photon!





Photosensors detect and transform the light produced by the scintillator into an **electrical signal**. This signal is proportional to the energy released inside the crystal by the interacting particle

Electronics & Analysis Software

The output pulses are proportional to the energies of the incident radiation, the ADC is used combined to a Multichannel Analyzer (MCA) to generate energy distributions (spectra) of radioactive samples




- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')**
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Potassium
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition
- 8) Passive Radon Measurements





Scintillator based on Lutetium (Lu) like LYSO (Lu_{1.8}Y_{0.2}SiO₅:Ce) has an high stopping power (high Z), high light yield and very short decay time (very fast signals).

2,6% of the natural Lu is ¹⁷⁶Lu, a radioisototope with a long half life decaying via two different beta decays followed by gamma emissions.

- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')
- 3) Background measurement (time run = 30')**
- 4) Calibration verification and tuning with Potassium Chloride sample (time run = 30')
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition
- 8) Passive Radon Measurements




- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Potassium Chloride sample (time run = 30')**
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition
- 8) Passive Radon Measurements




Fertilizer - Potassium Chloride

For our hands on take a 30 min spectra with the Potassium Chloride sample

Fertilizer




Lead shielding to reduce environmental background



U-238

- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Potassium
- 5) **Rock sample Spectrum (time run = 30')**
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition
- 8) Passive Radon Measurements

		Th-232	
99.98%	0.02%		
Pb214	26.8 minutes		53.2 (1.1%) 242.0 (7.46%) 295.2 (19.2%) 351.9 (37.1%) 785.9 (1.09%)
		At218	2 seconds
Bi214	19.7 minutes		609.3 (46.1%) 768.4 (4.89%) 806.2 (1.23%) 934.1 (3.16%) 1120.3 (15.0%) 1238.1 (5.92%) 1377.7 (4.02%) 1408.0 (2.48%) 1509.2 (2.19%) 1764.5 (15.9%)
Bi212	60.6 minutes		39.9 (1.1%) 727.3 (6.65%)
64.06%	35.94%		
Po212	304 nsec		
Tl208	3.1 minutes		277.4 (6.31%) 510.77 (22.6%) 583.2 (84.5%) 763.1 (1.81%) 860.6 (12.4%)



- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Potassium Chloride sample (time run = 30')
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification**
- 7) Analysis of spectra and superposition
- 8) Passive Radon Measurements



**Thorium
Lantern Mantle**



**Rare earth
uranium oxide**

Uranium Glazed Pottery




Uranium glass beads



**Decades-Old
Lenses**


Superimposed spectra



- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on LYSO crystal (time base = 10')
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Pot
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition**
- 8) Passive Radon Measurements

- 1) How to use the i-Spector Digital
- 2) Energy calibration of the system based on
- 3) Background measurement (time run = 30')
- 4) Calibration verification and tuning with Pot
- 5) Rock sample Spectrum (time run = 30')
- 6) Test sample radiation identification
- 7) Analysis of spectra and superposition
- 8) **Passive Radon Measurements (time run=60')**

Measurement of the adsorbed radon activity




Canisters of Activated Carbon

The amount of adsorbed radon by the canisters can be evaluated via the detection of the gamma rays emitted by the ^{214}Pb e dal ^{214}Bi .

Among the many available gamma emissions, the following nuclei are used as they are formed in a short time from the decay of Radon:

- ❖ 295 keV and 352 keV from ^{214}Pb
- ❖ 609 keV from ^{214}Bi




- SiPM is a High density (up to $10^4/\text{mm}^2$) matrix of diodes with a common output, working in Geiger-Müller regime
- Common bias is applied to all cells (few % over breakdown voltage)
- Each cell has its own quenching resistor (from 100kΩ to several MΩ)
- When a cell is fired an avalanche starts with a multiplicative factor of about 10^5 - 10^6
- The output is a fast signal ($t_{\text{rise}} \sim \text{ns}$; $t_{\text{fall}} \sim 50 \text{ ns}$) sum of signals produced by individual cells
- SiPM works as an analog photon detector


$$A_i \sim C (V_{\text{breakdown}} - V_{\text{bias}})$$

$$A = \sum A_i$$

- SiPM may be seen as a collection of binary cells, fired when a photon is absorbed
- “counting” cells provides an information about the intensity of the incoming light




The high uniformity of pixel structure guarantees no avalanche fluctuations



Linear response if the average number of photoelectrons/pixel is less than one

Number of pixel determines the SiPM **dynamic range**





Excellent Resolution


The **Dark Counts (DCR)** measure the rate at which a Geiger avalanche is randomly initiated by thermal emission.

Decrease DCR:

- lowering temperature
- lowering active volume – decrease V_{bias}
– small area



An avalanche generation can fire another cell by a photon; measuring the DCR for different thresholds is possible to define and evaluate the **Optical Cross Talk (OCT)** as:




Decrease OCT:

- low gain/ V_{bias}
- big pixel size
- trench for optical isolation

$$X_{talk} = \frac{DCR(1.5 ph)}{DCR(0.5 ph)}$$


After Pulse: It is a delayed avalanches triggered by the release of a charge carrier that has been produced in the original avalanche and has been trapped on an impurity



- Limiting the photon counting resolution
- The release of the trapped carriers is characterized by a typical decay time ~ 200 ns

The pulse amplitude depends on the pixel recovery state

$$\xi(\Delta t) = 1 - \exp(-\Delta t / \tau_r)$$



Decrease AP:

- low V_{bias}
- small pixel

Granite (intrusive igneous rock)

~1000 Bq/kg

**Tuff (magmatic rock)**

~1800 Bq/kg

**Marble (limestone rock)**


~100 Bq/kg

**Gypsum (sedimentary rock)**


~100 Bq/kg



Electromagnetic Radiation Spectrum



The human eye can see light radiation
from ~370 nm and ~700 nm



Il MeV: una grandezza conveniente per misurare le energie dei decadimenti radioattivi


Comunemente viene usato il Joule (J) come unità di misura per misurare l'energia. Durante i processi di decadimento l'energia rilasciata è trasportata dalle particelle è tipicamente dell'ordine di grandezza di 10^{-13} J.

E' comodo introdurre una unità di energia chiamata elettro-volt (eV) definita come l'energia guadagnata da un elettrone che attraversa una differenza di potenziale 1V.

$$1\text{eV} = 1.6 \times 10^{-19} \text{ C} \times 1 \text{ V} = 1.6 \times 10^{-19} \text{ Joule}$$

ossia, anche

$$1 \text{ Joule} = 0.6 \times 10^{19} \text{ eV}$$




Nei processi nucleari le energie in gioco, dell'ordine di 10^{-13} J, sono dunque dell'ordine di

$$1 \text{ MeV (mega elettro-volt)} = 10^6 \text{ eV}$$

Cosa significa che due misure sono compatibili?

Per confrontare due misure della stessa grandezza (in questo caso le abbondanze di U, Th e K) occorre confrontare gli intervalli di confidenza. Le due misure si dicono compatibili se i rispettivi intervalli di incertezza hanno intersezione non nulla.



$$aK[\%]_1 = 1.5 \pm 0.2 \text{ e } aK[\%]_2 = 1.8 \pm 0.2$$

$$aK[\%]_3 = 1.2 \pm 0.2 \text{ e } aK[\%]_4 = 1.7 \pm 0.2$$

Compatibili




Non compatibili



Luminescenza: Materiali luminescenti assorbono energia e la riemettono sotto forma di luce visibile.

Fluorescenza: l'emissione avviene subito dopo l'assorbimento (10^{-8} s).

Fosforescenza: l'emissione è ritardata (lo stato eccitato è metastabile). In questo caso il tempo fra l'assorbimento e la riemissione può durare dai μ s alle ore (dipende dal materiale).



In prima approssimazione l'evoluzione temporale del processo di riemissione può essere descritto da una semplice legge esponenziale:

$$N = \frac{N_0}{\tau_d} \cdot e^{-t/\tau_d}$$


Dove N è il numero di fotoni emessi al tempo t , τ_d il tempo di decadimento ed N_0 il numero totale di fotoni emessi.

Il tempo di assorbimento dell'energia (eccitazione degli atomi e delle molecole) è in generale molto più breve del tempo di riemissione.

Scintillatori Inorganici

Il meccanismo di scintillazione negli scintillatori inorganici è caratteristico della struttura a bande elettroniche che si trovano nei cristalli.

Quando una particella entra in un cristallo possono accadere 2 processi:



- i. si eccita un elettrone dalla banda di valenza in quella di conduzione, creando così un elettrone ed una lacuna liberi. (ionizzazione)
- ii. si crea un eccitone spostando un elettrone dalla banda di valenza in quella degli eccitoni (posta appena al di sotto della banda di conduzione). In questo caso elettrone e lacuna rimangono legati, ma possono muoversi liberamente (in coppia) nel cristallo. (eccitazione)


Se il cristallo contiene delle impurità (sono necessarie), si possono creare localmente dei livelli elettronici nella banda delle energie proibite. Gli atomi di impurità sono chiamati centri attivatori.

Se una lacuna libera od una lacuna di un eccitone incontra uno di questi centri attivatori, può ionizzare uno di questi atomi di impurità. Se ora arriva un altro elettrone, questo cade nel buco (lacuna) lasciato dalla ionizzazione precedente → si emette luce (se tale modo di diseccitazione è permesso).

Scintillatori Organici

Gli scintillatori organici sono dei composti di idrocarburi. In questi composti la luce di scintillazione deriva da transizioni degli elettroni di valenza liberi delle molecole.

L'eccitazione verso stati di energia superiore a S_1 è seguita da una rapidissima transizione (nell'ordine del ps) non radiativa verso S_1 . Da questo livello le molecole si diseccitano verso lo stato fondamentale S_0 o uno degli stati vibrazionali emettendo luce di fluorescenza.



Caratteristiche generali degli scintillatori

Cristalli alcalino-alogeno: NaI, CsI, BGO, BaF₂

Buona linearità, lenti tempi di risposta


Spettroscopia gamma




Scintillatori organici: antracene, stilbene

Buona velocità, bassa resa in luce

Spettroscopia beta e per neutroni veloci



→ Scintillatore




A cosa servono gli scintillatori?

- A misurare l'energia ceduta
- A misurare il tempo di passaggio delle radiazioni


CsI(Tl)

Lo Ioduro di Cesio drogato al Tallio, CsI(Tl) , è un materiale scintillatore che ha un'alta uscita in luce ed è relativamente lento. Il tempo di decadimento del cristallo, di circa 1100 ns.

Densità	4.51 g/cm ³
Indice di rifrazione per λ max	1.788
Punto di ebollizione	620 °C
Uscita in luce	95%
Costante di decadimento	1100 ns
λ del max di emissione	580 nm
Principali applicazioni	rivelazione γ , particelle pesanti



Dai fotoni alla corrente elettrica



La luce viene convertita in segnale elettrico mediante l'impiego di un fotomoltiplicatore

L'intensità dell'impulso di corrente in uscita è proporzionale all'energia del fotone incidente

Il canale di acquisizione è proporzionale all'energia dei fotoni incidenti