


B.1.9

SG6119B

Activity of the ^{60}Co 

Dedicated kit	
Description	pp.
SP5600EMU Emulation Kit	180



Requirements

Gamma Radioactive Source 

Difficulty



Execution Time



Data Analysis

NO

Radioactive Sources

YES

Equipment

SP5600EMU - Emulation Kit

Model	DT4800	DT5770
Description	Digital Detector Emulator	Digital Multi-Channel Analyzer
		
	p. 194	p. 194

Purpose of the experiment

Determine the activity of a ^{60}Co source from its gamma spectrum.
Learn about the meaning of the sum peak, visible in the spectrum of some radioactive sources.

See the Application



Fundamentals

The ^{60}Co spectrum presents two distinct gamma photopeak in its spectrum, respectively corresponding to photons γ_1 and γ_2 at 1.17 MeV and 1.33 MeV. For the purpose of this experiment, we can assume that each of these gamma rays are isotropically distributed. In other words, if γ_1 departs in a particular direction, γ_2 can go in any direction that it wishes. There is a certain probability that γ_2 will go in the same direction as γ_1 . If this occurs the energies of γ_1 and γ_2 will be summed in the detector. Hence a sum peak will show up in the spectrum, at nearly 2.5 MeV.

We can estimate the activity of the source by calculating the counts under the two main peaks and under the sum peak, i.e. calculating their area Σ . For the case of ^{60}Co , we have that the counts under the sum peak can be evaluated as


$$\Sigma(\text{SUM}) = \frac{\Sigma_1 \Sigma_2}{A \tau}$$

Where A is the activity of the source and t is the acquisition time.

Therefore, fitting the peaks with a gaussian and calculating their area, it is possible to estimate the activity of the ^{60}Co source used to record the available spectrum.

Carrying out the experiment

To perform the experiment, connect the DT4800 output to the input channel of the MCA DT5770 and use the DT4800 GPO as digitizer “trigger IN”. The Emulation Control Software Interface allows user to generate exponential decay signals with programmable rise time and fall time and it is possible to emulate signals from ^{137}Cs radioactive. The spectrum can be recorded and analyzed by the MCA.



Experimental setup block diagram for the experiment.

Results

The student should verify that, after the spectrum calibration, the sum peak is nearly at 2.5 MeV. From the formula given above, using the live time in seconds, the student can estimate the activity of ^{60}Co directly in Bq. A calculation made for a spectrum acquired in 100 seconds gives an activity of nearly 264 kBq.



The ^{60}Co complete spectrum acquired by the MCA DT5770 and plotted by the Emulation Software

This experiment is also possible with the following kits

A SP5600C Educational kits



see
p. 179

B SP5600AN Educational kits ★ PREMIUM



see
p. 179

C SP5640 GammaEDU



see
p. 183

D SP5630EN Environmental kit



see
p. 181

E SP5600ENP Environmental kit Plus



see
p. 181