

C.1.5 Random Coincidence

SG6213D



Dedicated kit	
Description	pp.
SP5620CH Cosmic Hunter	182



Requirements

No other tools are needed.

Difficulty	Execution Time	Data Analysis	Radioactive Sources
    	   	YES	YES

Equipment

SP5620CH - Cosmic Hunter

Model	SP5621	SP5622 (x2)
Description	Coincidence Module  p. 182	Detection System  p. 187

Purpose of the experiment

To understand the potential for accidental counts coming from double tile coincidence.

See the
Application



Fundamentals

Once the geometry of the detectors has been defined and coincidence is confirmed it becomes important to estimate the number of random coincidences. Random coincidences derive from simultaneous or nearly-simultaneous pulses caused by accidental discharges (i.e. noise), and not by particles with a trajectory within the volume determined by the solid angle of the geometry. The probability that a particle crosses the detector is a function of the surface of the tile itself and of the average rate. Therefore, the probability of a random coincidence is proportional to the pulse duration.


An evaluation of random coincidence contribution [R_{random}] can be obtained by a simple theoretical calculation, Janossy method based Error! Reference source not found:

$$R_{random} = 2 \cdot R_A \cdot R_C \cdot \tau$$

Where τ is the pulse duration (700 ns), and R_A and R_C are the event rate of each scintillating tile.

Carrying out the experiment

Connect the cable connectors of the two SP5622 to the tile inputs located on the rear panel of the SP5621 module. Power on the SP5621 module and start the acquisition via the front panel START button. When a charged particle crosses the black tile it's energy is converted into scintillation light. The photons which are produced are detected by the photosensor and converted into an electrical signal. The number of counts for each scintillator may be viewed via the SP5621 display. Select the scintillators coincidence via the related button on the front panel, then select the integration time of the measurement.




Experimental setup block diagram.

Before starting acquisition choose the system geometry. Be sure to keep this geometry constant for the duration of the experiment. Take and record more data to obtain statistical significance.

Results

Double tile coincidence plays a key role in a great many Physics experiments. The random coincidence rate allows you to evaluate the data quality via estimation of the Signal to Background ratio [SBR].



Trend of the Count Rate and Random Rate as a function of the time. The plot on the left side is an enlargement of the main plot and underlines the deviation between the measured coincidence rate and the real one, obtained via the random rate subtraction.

This experiment is also possible with the following kits

A SP5600D Educational kits



see
pp. 182,
179

B SP5600AN Educational kits ★ PREMIUM



see
p. 179